

Safety Data Sheet

Issue Date: 23-Aug-2013

Revision Date: 08-May-2018

Version 2

1. IDENTIFICATION

Product Identifier

Product Name ProTech Prosthetic Cleaner

Other means of identification

SDS # PROTECH-001

Recommended use of the chemical and restrictions on use

Recommended Use Prosthetic cleaner.

Details of the supplier of the safety data sheet

Supplier Address

ProTech Professional Products, Inc.
2900 N.W. Commerce Park Dr.
Suite 10
Boynton Beach, FL 33426

Emergency Telephone Number

Company Phone Number 1-800-872-8898
1-561-493-9818
Emergency Telephone (24 hr) INFOTRAC 1-352-323-3500 (International)
1-800-535-5053 (North America)

2. HAZARDS IDENTIFICATION

Appearance White, free-flowing solid

Physical state Solid

Odor Odorless

Classification

Acute toxicity - Oral	Category 4
Skin corrosion/irritation	Category 2
Serious eye damage/eye irritation	Category 2

Signal Word

Warning

Hazard statements

Harmful if swallowed
Causes skin irritation
Causes serious eye irritation



Precautionary Statements - Prevention

Wash face, hands and any exposed skin thoroughly after handling
Do not eat, drink or smoke when using this product
Wear protective gloves/protective clothing/eye protection/face protection

Precautionary Statements - Response

IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing
 If eye irritation persists: Get medical advice/attention
 IF ON SKIN: Wash with plenty of water
 Take off contaminated clothing and wash it before reuse
 If skin irritation occurs: Get medical advice/attention
 IF SWALLOWED: Call a POISON CENTER or doctor/physician if you feel unwell
 Rinse mouth

Precautionary Statements - Disposal

Dispose of contents/container to an approved waste disposal plant

Other hazards

Harmful to aquatic life with long lasting effects

3. COMPOSITION/INFORMATION ON INGREDIENTS

Chemical Name	CAS No.	Weight-%
Trade Secret	Proprietary	80-90

The specific chemical identity of this composition is being withheld as a trade secret

4. FIRST AID MEASURES

First Aid Measures

- Eye Contact** Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Seek immediate medical attention/advice.
- Skin Contact** Wash with soap and water. Take off contaminated clothing. Wash contaminated clothing before reuse. If skin irritation persists, call a physician.
- Inhalation** Remove to fresh air. If not breathing, give artificial respiration. Get medical attention immediately.
- Ingestion** Rinse mouth. Do NOT induce vomiting. Give water to conscious/alert person. Never give anything by mouth to an unconscious person. Get medical attention immediately.

Most important symptoms and effects

Symptoms Irritation to lungs and respiratory tract if dust is inhaled. Eye contact may result in redness, pain, blurred vision, burning sensation. Prolonged contact may even cause severe skin irritation or mild burn. Nausea. Vomiting.

Indication of any immediate medical attention and special treatment needed

Notes to Physician Treat symptomatically.

5. FIRE-FIGHTING MEASURES

Suitable Extinguishing Media

Use extinguishing measures that are appropriate to local circumstances and the surrounding environment.

Unsuitable Extinguishing Media Not determined.

Specific Hazards Arising from the Chemical

Storage of large masses of this product can trap heat and lead to ignition of paper bags. Grinding or intensive mixing can cause ignition of oxidizable material present.

Protective equipment and precautions for firefighters

As in any fire, wear self-contained breathing apparatus pressure-demand, MSHA/NIOSH (approved or equivalent) and full protective gear.

6. ACCIDENTAL RELEASE MEASURES

Personal precautions, protective equipment and emergency procedures

Personal Precautions Wear protective equipment including rubber boots, rubber gloves, rubber apron, and chemical goggles.

Environmental precautions

Environmental precautions Prevent material from entering surface waters, drains or sewers and open soil.

Methods and material for containment and cleaning up

Methods for Containment Prevent further leakage or spillage if safe to do so.

Methods for Clean-Up Sweep up and shovel into suitable containers for disposal. Clean up in accordance with all applicable regulations. Discard any product, residue, disposable container or liner in full compliance with federal, state, and local regulations.

7. HANDLING AND STORAGE

Precautions for safe handling

Advice on Safe Handling Wash thoroughly after handling. Do not eat, drink or smoke when using this product. Use personal protection recommended in Section 8. Take any precaution to avoid mixing with combustibles. Avoid contact with skin, eyes or clothing. Emptied container retains product residue. Observe all labeled safeguards until container is cleaned, reconditioned or destroyed.

Conditions for safe storage, including any incompatibilities

Storage Conditions Keep containers tightly closed in a dry, cool and well-ventilated place. Keep bags dry at all times. Closely stacked bags should not exceed a 4 foot cube. Do not store near combustible materials.

Incompatible Materials Compounds containing halides or active halogens can release the respective halogen if moisture is present; Mixture with cyanides can release poisonous hydrogen gas; Mixture with heavy metal salts can release flammable oxygen gas.

8. EXPOSURE CONTROLS/PERSONAL PROTECTION

Exposure Guidelines No exposure limits noted for ingredient(s).

Appropriate engineering controls

Engineering Controls Ventilation systems. Eyewash stations. Showers.

Individual protection measures, such as personal protective equipment

Eye/Face Protection Chemical goggles or equivalent eye protection.

Skin and Body Protection Leather gloves. Normal work clothing (long sleeved shirts and long pants) is recommended.

Respiratory Protection NIOSH-approved dust respirator or mask in the absence of adequate environmental controls at the point of use.

General Hygiene Considerations Handle in accordance with good industrial hygiene and safety practice.

9. PHYSICAL AND CHEMICAL PROPERTIES

Information on basic physical and chemical properties

Physical state	Solid	Odor	Odorless
Appearance	White, free-flowing solid	Odor Threshold	Not determined
Color	White		

<u>Property</u>	<u>Values</u>	<u>Remarks • Method</u>
pH	Not determined	
Melting point / freezing point	Not determined	
Boiling Point / Boiling Range	Decomposes	
Flash Point	None	Tag Closed Cup
Evaporation Rate	+1	(butyl acetate = 1)
Flammability (Solid, Gas)	Not determined	
Flammability Limit in Air		
Upper Flammability Limit	Not applicable	
Lower Flammability Limit	Not applicable	
Vapor Pressure	Nil	
Vapor Density	Not applicable	
Relative Density	1.1-1.4	(1=Water)
Water Solubility	25.6%	
Solubility in other solvents	Not determined	
Partition Coefficient	Not determined	
Autoignition temperature	Not applicable	
Decomposition Temperature	Not determined	
Kinematic Viscosity	Not determined	
Dynamic Viscosity	Not determined	
Explosive Properties	Not determined	
Oxidizing Properties	Not determined	

10. STABILITY AND REACTIVITY

Reactivity

Not reactive under normal conditions.

Chemical Stability

Stable under recommended storage conditions.

Possibility of Hazardous Reactions

See below - Incompatible Materials.

Hazardous Polymerization Hazardous polymerization does not occur.

Conditions to Avoid

None known.

Incompatible Materials

Compounds containing halides or active halogens can release the respective halogen if moisture is present; Mixture with cyanides can release poisonous hydrogen gas; Mixture with heavy metal salts can release flammable oxygen gas.

Hazardous Decomposition Products

Oxygen.

11. TOXICOLOGICAL INFORMATION

Information on likely routes of exposure

Product Information

- Eye Contact** Causes serious eye irritation.
- Skin Contact** Causes skin irritation. May be harmful in contact with skin.
- Inhalation** Avoid inhalation of dust.
- Ingestion** Harmful if swallowed.

Component Information

Chemical Name	Oral LD50	Dermal LD50	Inhalation LC50
Trade Secret	= 1204 mg/kg (Rat)	> 11000 mg/kg (Rabbit)	> 14 mg/L (Rat) 1 h

Information on physical, chemical and toxicological effects

- Symptoms** Please see section 4 of this SDS for symptoms.

Delayed and immediate effects as well as chronic effects from short and long-term exposure

- Carcinogenicity** Based on the information provided, this product does not contain any carcinogens or potential carcinogens as listed by OSHA, IARC or NTP.

Numerical measures of toxicity

The following values are calculated based on chapter 3.1 of the GHS document .

- ATEmix (oral)** 1,480.00 mg/kg
- ATEmix (dermal)** 2,459.00 mg/kg
- ATEmix (inhalation-vapor)** 6.00 mg/L

12. ECOLOGICAL INFORMATION

Ecotoxicity

Harmful to aquatic life. Harmful to aquatic life with long lasting effects.

Persistence/Degradability

Not determined.

Bioaccumulation

Not determined.

Mobility

Not determined

Other Adverse Effects

Not determined

13. DISPOSAL CONSIDERATIONS

Waste Treatment Methods

- Disposal of Wastes** Disposal should be in accordance with applicable regional, national and local laws and regulations.
- Contaminated Packaging** Disposal should be in accordance with applicable regional, national and local laws and regulations.

14. TRANSPORT INFORMATION

Note Please see current shipping paper for most up to date shipping information, including exemptions and special circumstances.

DOT Not regulated

IATA Not regulated

IMDG Not regulated

15. REGULATORY INFORMATION

International Inventories

Chemical Name	TSCA	DSL/NDSL	EINECS/E LINC S	ENCS	IECSC	KECL	PICCS	AICS
Trade Secret	X	X	X		X	X	X	X

Legend:

- TSCA - United States Toxic Substances Control Act Section 8(b) Inventory*
- DSL/NDSL - Canadian Domestic Substances List/Non-Domestic Substances List*
- EINECS/ELINCS - European Inventory of Existing Chemical Substances/European List of Notified Chemical Substances*
- ENCS - Japan Existing and New Chemical Substances*
- IECSC - China Inventory of Existing Chemical Substances*
- KECL - Korean Existing and Evaluated Chemical Substances*
- PICCS - Philippines Inventory of Chemicals and Chemical Substances*
- AICS - Australian Inventory of Chemical Substances*

US Federal Regulations

SARA 313

Not determined

US State Regulations

U.S. State Right-to-Know Regulations

Not determined

16. OTHER INFORMATION

<u>NFPA</u>	Health Hazards	Flammability	Instability	Special Hazards
	Not determined	Not determined	Not determined	Not determined
<u>HMIS</u>	Health Hazards	Flammability	Physical hazards	Personal Protection
	Not determined	Not determined	Not determined	Not determined

Issue Date: 23-Aug-2013
Revision Date: 08-May-2018
Revision Note: Regulatory update

Disclaimer

The information provided in this Safety Data Sheet is correct to the best of our knowledge, information and belief at the date of its publication. The information given is designed only as a guidance for safe handling, use, processing, storage, transportation, disposal and release and is not to be considered a warranty or quality specification. The information relates only to the specific material designated and may not be valid for such material used in combination with any other materials or in any process, unless specified in the text.

End of Safety Data Sheet